

# Detection of *Candidatus Liberibacter solanacearum* in Carrot Seed

**DECEMBER 2023**

**Developed by ISHI**

All rights reserved - ©2023 ISF

## Detection of *Candidatus Liberibacter solanacearum* in Carrot Seed

|                     |   |
|---------------------|---|
| <b>Crop:</b>        | <i>Daucus carota</i> (carrot)               |
| <b>Pathogen(s):</b> | <i>Candidatus Liberibacter solanacearum</i> |
| <b>Version:</b>     | 2.2 (December 2023)                         |

### BACKGROUND

In ISHI's view, there is a lack of evidence for seed transmission of *Candidatus Liberibacter solanacearum* (*Ca. L. solanacearum*) via carrot seed. Seed transmission has been reported in a single publication (Bertolini *et al.*, 2015). Several other studies have challenged this finding (Haapalainen, 2014; Loiseau *et al.*, 2017a and 2017b; Oishi *et al.*, 2017; Mawassi *et al.*, 2018; Fujikawa *et al.*, 2020).

As a consequence of the report by Bertolini *et al.* (2015), several governments started to impose emergency measures and import restrictions on seed of carrot and other *Apiaceae* species. Lately, however, New Zealand's national plant protection organisation reported that its full technical assessment found no new scientific evidence to support the finding by Bertolini *et al.* (2015), that *Ca. L. solanacearum* is seed transmitted on the *Apiaceae* family. As a result, it no longer considers *Apiaceae* (i.e., carrot, celery, chervil, fennel, parsley and parsnip) seed to be a risk pathway for the introduction of *Ca. L. solanacearum*. As a trade-enabling measure New Zealand removed the *Apiaceae* specific requirements from its Import Health Standard (see <https://www.mpi.govt.nz/dmsdocument/1151/direct>).

To facilitate trade of carrot seed, ISHI has developed a detection method based on real-time PCR, a requirement from some governments to obtain a Phytosanitary Certificate for exporting seed. ISHI calls on researchers seeking to develop detection methods for *Ca. L. solanacearum* to use the primer sets used in this method to harmonise methods used for phytosanitary certification. These primers and probes have been developed by Li *et al.* (2009) and are also recommended in EPPO's Diagnostic Standard PM 7/143 (1) for the detection of *Ca. L. solanacearum* in both asymptomatic and symptomatic plant material or in vectors (see <https://onlinelibrary.wiley.com/doi/epdf/10.1111/epp.12611>).

According to ISPM 38, *molecular and serological diagnostic methods are considered indirect protocols to detect pests in seeds. These methods may give a positive result even when no viable pests are present. Consequently, when testing seeds with these methods, results should be interpreted carefully. Confirmatory tests or additional tests based on a different biological principle may be required to confirm the presence of a viable pest in a sample* (see <https://www.ippc.int/en/publications/84340/>).

A qPCR assay, an indirect test, detects both infectious and non-infectious bacterial particles as well as DNA from viable and non-viable bacteria (see [http://www.worldseed.org/wp-content/uploads/2015/10/Indirect\\_Seed\\_Health\\_Tests\\_2013.pdf](http://www.worldseed.org/wp-content/uploads/2015/10/Indirect_Seed_Health_Tests_2013.pdf)). A positive SE-qPCR is, therefore, only indicative of a suspected lot. A confirmatory test that determines the presence of


viable and pathogenic *Ca. L. solanacearum* is needed to confirm the health status of the lot. Such a confirmation test, usually a pathogenicity or bioassay, is currently not available for *Ca. L. solanacearum* due to the uncultivable nature of this pathogen. Furthermore, as efforts to show seed transmission have not been successful, it is not possible to develop a direct test that confirms viability of the pathogen and its pathogenicity.

The qPCR has been designed to identify seed lots that are *not* infected with the target pathogen and validated to show there are no false negative results obtained.

## PRINCIPLE

Detection of *Ca. L. solanacearum* on carrot seed is assessed by a seed extract qPCR assay (SE-qPCR). Extraction of bacteria located in the phloem sieve tubes of the carrot seed coat is enhanced by stomaching. The test is complete if no *Ca. L. solanacearum* is detected and the seed lot is considered healthy. A positive SE-qPCR is only indicative of the presence of DNA sequence of *Ca. L. solanacearum* and gives no information on viability and pathogenicity.

The method process workflow is presented in Figure 1.



**Figure 1.** Method process workflow.

## METHOD VALIDATION

This method has been validated by ISHI (Chitrampalam *et al.*, 2020). The validation included an inter-laboratory comparative test. The performance of the two main components of the ISHI assay, namely the Qiagen DNA extraction and the Li *et al.* (2009) real-time PCR method, has also been validated independently of ISHI for analytical specificity, sensitivity and repeatability (Loiseau *et al.*, 2015).

## RESTRICTIONS ON USE

Before using this protocol routinely, it is necessary to verify its performance, especially when material and consumables from different suppliers are used. Technical details on the reagents/material used in the validation study (e.g., supplier's information) are provided in the protocol and the validation report.

This test method is suitable for untreated seed.

This test method has not been validated for seed treated with protective chemicals or biological substances. If a user chooses to test treated seed using this method, it is the responsibility of the user to determine empirically (through analysis, sample spiking, or experimental comparisons) whether the protective chemicals or biological substances have an effect on the method results.

## METHOD EXECUTION

To ensure process standardization and valid results, it is strongly recommended to follow the [Best Practices for Seed Health Tests](#) developed by ISHI.

## SAMPLE AND SUBSAMPLE SIZE

The recommended minimum sample size is 20,000 seeds with a maximum subsample size of 10,000 seeds.

NOTE: ISHI is not aware of any epidemiological information that warrants seed samples of any particular size.

## REVISION HISTORY

| Version | Date          | Changes (minor editorial changes not indicated)   |
|---------|---------------|---|
| 1.0     | January 2015  | First version of the protocol.  |
| 1.1     | July 2016     | Additional information from new publications added, including references to the new publications.<br>The use of PMA has been deleted as its efficacy is less than 100%. Positive PCR results can be obtained due to incomplete removal of DNA from dead bacteria. |
| 2.0     | November 2019 | An internal amplification control is added to the SE-qPCR.  |
| 2.1     | June 2020     | Background section updated to latest information.   |
| 2.2     | December 2023 | Protocol updated in accordance with ISHI protocol guidelines.   |

## Protocol for detecting *Candidatus Liberibacter solanacearum* in Carrot Seed

### I. DETECTION OF TARGET DNA BY SEED-EXTRACT qPCR

For PCR methods, in-house method optimization is often necessary, see [Best Practices for PCR Assays in Seed Health Tests](#) for information.

#### Materials

- Stomacher and stomacher bags
- DNA isolation kit e.g., Qiagen DNeasy Plant Mini kit
- qPCR mix, primers (Table I.1) and PCR equipment
- Controls (Table I.2)
- PBS buffer (Table I.3)
- Centrifuge

**Table I.1.** Primer sequences and references.

| Name     | Sequence (5' – 3')                         | Source                         |
|----------|--|--------------------------------|
| Lso-F    | GTC GAG CGC TTA TTT TTA ATA GGA'           | Li <i>et al.</i> , 2006, 2009  |
| HLBr     | GCG TTA TCC CGT AGA AAA AGG TAG            |                                |
| HLBp     | FAM – AGA CGG GTG AGT AAC GCG – BHQ1       |                                |
| Acat2-F  | TGT AGC GAT CCT TCA CAA G                  | Koenraadt <i>et al.</i> , 2014 |
| Acat2-R  | TGT CGA TAG ATG CTC ACA AT                 |                                |
| Acat1-Pr | HEX – CTT GCT CTG CTT CTC TAT CAC G – BHQ2 |                                |

**Table I.2.** Types of controls used.

| Control type                            | Description   |
|---|---|
| Positive extraction control (PEC)       | <i>Acidovorax cattleyae</i> spike   |
| Negative extraction control (NEC)       | PBS buffer  |
| Internal amplification control (IAC)    | <i>Acidovorax cattleyae</i> spike   |
| Positive process control (PPC)          | Seeds infected with <i>Ca. L. solanacearum</i> , or healthy seed lots spiked with <i>Ca. L. solanacearum</i> extracts |
| Negative process control (NPC)          | Healthy seeds ( <i>Ca. L. solanacearum</i> negative seeds)  |
| Positive amplification control 1 (PAC1) | <i>Ca. L. solanacearum</i> DNA  |
| Positive amplification control 2 (PAC2) | <i>Acidovorax cattleyae</i> DNA   |
| Non template control (NTC)              | Nucleic acid free water   |

**Table I.3.** Phosphate buffered saline (PBS) buffer.

| Compound                         | Amount/L |
|----------------------------------|----------|
| NaCl                             | 8.00 g   |
| KH <sub>2</sub> PO <sub>4</sub>  | 0.24 g   |
| Na <sub>2</sub> HPO <sub>4</sub> | 1.44 g   |
| KCl                              | 0.20 g   |

## 1. Extraction of bacteria from seeds

- 1.1. Add the subsamples, the positive process control (PPC) and the negative process control (NPC) to sterile stomacher bags with 100 mL sterile PBS buffer (ratio of 100 mL buffer/10,000 seeds).
- 1.2. As negative extraction control (NEC), include a sterile stomacher bag with 100 mL sterile PBS buffer only (without seeds).
- 1.3. As positive extraction control (PEC) and internal amplification control (IAC), add 100 µL of Acat spike (OD<sub>600</sub>=0.6) to each of the subsamples, the PPC, the NPC and the NEC.
- 1.4. Macerate all samples and controls (including NEC without seeds) in a stomacher machine until all the seeds are crushed.

Note: With InterScience Bag mixer good results were obtained by stomaching for 4 min at speed 3 and gap 5.

- 1.5. Pipet 2 mL of the extract from the filtered side of the stomacher bag into a 2 mL tube.

## 2. DNA extraction

- 2.1. Centrifuge the 2 mL filtered extract to pellet large seed debris (500 × *g*, 1 min) and transfer 1 mL of the supernatant to a new 1.5 mL micro-centrifuge tube.
- 2.2. Centrifuge at 10,000 × *g* for 5 min to pellet the bacteria. Remove the supernatant and continue the DNA extraction with the Qiagen DNeasy kit on the pellet in each tube.
- 2.3. Check that the salts in the Qiagen DNeasy Plant Mini Kit are fully suspended in the AP1 buffer, if necessary, place at 65 °C, until the salts are dissolved.
- 2.4. Add 400 µL Buffer AP1 and 4 µL RNase-A stock solution (100 mg/mL) to the pellet. Re-suspend the pellet, vortex and incubate at 65 °C for 10 min.
- 2.5. Add 130 µL Buffer P3 to the lysate, mix, and incubate for 5 min on ice or at +4 °C.
- 2.6. Centrifuge the lysate at 20,000 × *g* for 5 min.
- 2.7. Pipet the lysate into the QIA shredder Mini spin column (lilac) placed in a 2 mL collection tube, and centrifuge at 20,000 × *g* for 2 min.
- 2.8. Transfer the flow-through into a new 1.5 mL tube without disturbing the cell-debris pellet.
- 2.9. Add 1.5 volumes of Buffer AW1 to the cleared lysate and mix by pipetting.

- 2.10. Pipet 650  $\mu\text{L}$  of the mixture, including any precipitate that may have formed, into the DNeasy Mini spin column (white) placed in a 2 mL collection tube and centrifuge at  $6,000 \times g$  for 1 min and discard the flow-through.
- 2.11. Add 500  $\mu\text{L}$  Buffer AW2 and centrifuge at  $6,000 \times g$  for 1 min and discard the flow-through.
- 2.12. Repeat the wash step by adding 500  $\mu\text{L}$  Buffer AW2 to the DNeasy Mini spin column, and centrifuge at  $20,000 \times g$  for 2 min to dry the membrane.
- 2.13. Elute the extracted and purified DNA by adding 100  $\mu\text{L}$  AE Buffer, incubate at room temperature for 5 min and centrifuge into a clean 1.5 mL tube at  $6,000 \times g$  for 1 min.

Note: The assay has been validated with the Qiagen DNeasy Plant kit. If a different DNA isolation kit is used, it is necessary to verify its performance.

### 3. qPCR

Note: Work on ice whenever possible and minimize probes exposure to light.

- 3.1. Test the subsample extracts, NTC, PPC, NPC, NEC, PAC1 and PAC2 with the duplex Lso-Acat qPCR. Table I.4 provides an example for the reaction mixtures. However, reaction mixture and conditions need to be checked and/or optimized within each laboratory.
- 3.2. Add *in duplicate* the subsample extracts, NTC, PPC, NPC, NEC, PAC1 and PAC2 to the PCR reaction mixture in a suitable PCR reaction tube or plate.
- 3.3. Run the assay based on the conditions described in Table I.5 for the duplex qPCR.

Note: If different qPCR mixtures and amplification programs are used, it is necessary to verify their performance.

**Table I.4.** Example qPCR mix for the duplex Lso-Acat qPCR.

| Component                                  | Per reaction (in $\mu\text{L}$ ) | Final concentration |
|--|----------------------------------|---------------------|
| PCR grade $\text{H}_2\text{O}$             | 4.51                             |                     |
| qPCR Master Mix (2 $\times$ ) <sup>a</sup> | 7.50                             | 1 $\times$          |
| Lso-F (25 $\mu\text{M}$ )                  | 0.18                             | 0.3 $\mu\text{M}$   |
| HLBr (25 $\mu\text{M}$ )                   | 0.18                             | 0.3 $\mu\text{M}$   |
| HLBp (25 $\mu\text{M}$ )                   | 0.12                             | 0.2 $\mu\text{M}$   |
| Acat2-F (25 $\mu\text{M}$ )                | 0.18                             | 0.3 $\mu\text{M}$   |
| Acat2-R (25 $\mu\text{M}$ )                | 0.18                             | 0.3 $\mu\text{M}$   |
| Acat1-Pr (10 $\mu\text{M}$ )               | 0.15                             | 0.1 $\mu\text{M}$   |
| Sample                                     | 2.00                             |                     |
| <b>Total</b>                               | <b>15.00</b>                     |                     |

<sup>a</sup> TaqMan™ Universal Master Mix II (Applied Biosystems) for example.

**Table I.5.** PCR conditions for the duplex Lso-Acat qPCR.

| Step      | Temperature | Duration |
|-----------|-------------|----------|
| hold      | 95 °C       | 10 min   |
| 40 cycles | 95 °C       | 15 sec   |
|           | 60 °C       | 60 sec   |

#### 4. Interpretation and decisions

Cut-off values must be established by each laboratory for their positive and internal amplification controls prior to the assay being used on routine samples. For recommendations on setting cut-off values, see [Real-time PCR, an 'indirect' test used for pre-screening in seed health methods](#).

The cut-off for the IAC *Acidovorax cattleyae* (Acat) is set preferably below Cq 35. A cut-off value of Cq 40 for the Lso qPCR was used for the validation study. See Table I.6 for the interpretation of the qPCR results.

Test results are only valid when all included controls presented in Table I.2 give the expected results.

**Table I.6.** Interpretation and decision table duplex Lso-Acat qPCR.

| Lso qPCR | Acat qPCR (IAC) | qPCR result   |
|----------|-----------------|---|
| Cq < 40  | Not relevant    | Target DNA for <i>Ca. L. solanacearum</i> detected    |
| Cq ≥ 40  | Cq ≥ 35         | Inconclusive, repeat PCR                              |
| Cq ≥ 40  | Cq < 35         | No target DNA for <i>Ca. L. solanacearum</i> detected |

## REFERENCES

- Bertolini, E., Teresani, G. R., Loiseau, M., Tanaka, F. A. O., Barbé, S., Martínez, C., Gentit, P., López, M. M. and Cambra, M. (2015). Transmission of '*Candidatus Liberibacter solanacearum*' in carrot seeds. *Plant Pathology*, **64**(2), 276–285.
- Chitrampalam, P., Asma, M. and Hiddink, G. (2020). Detection of *Candidatus Liberibacter solanacearum* on Carrot Seed by Seed Extract PCR. Validation report, International Seed Federation (ISF), Nyon, Switzerland. <https://worldseed.org/our-work/seed-health/ishi-method-development-and-validation/>.
- EPPO (2020). PM 7/143 (1) '*Candidatus Liberibacter solanacearum*'. *EPPO Bulletin*, **50**(1), 49–68. <https://doi.org/10.1111/epp.12611>.
- Fujikawa, T., Yamamura, K., Osaki, K., Onozuka, N., Taguchi, M., Sasaki, A., and Sato, M. (2020). Seed transmission of '*Candidatus Liberibacter solanacearum*' is unlikely in carrot. *Journal of General Plant Pathology*, **86**, 266–273. <https://doi.org/10.1007/s10327-020-00927-1>.
- Haapalainen, M. (2014). Biology and epidemics of *Candidatus Liberibacter* species, psyllid-transmitted plant-pathogenic bacteria. *Annals of Applied Biology*, **165**, 172–198.



- Koenraad, H., van Vliet, A., Jodłowska, A., Woudt, B., Ebskamp, M. and Bruinsma, M. (2014). Detection of *Acidovorax citrulli* in seed extracts of melon and watermelon with TaqMan PCR. 7th ISTA Seed Health Symposium. Edinburgh, United Kingdom, 12–14 June 2014.
- Li, W., Abad, J. A., French-Monar, R. D., Rascoe, J., Wen, A., Gudmestad, N. C., Secor, G. A., Lee, I. M., Duan, Y., and Levy, L. (2009). Multiplex real-time PCR for detection, identification and quantification of ‘*Candidatus Liberibacter solanacearum*’ in potato plants with zebra chip. *Journal of Microbiological Methods*, **78**, 59–65.
- Li W., Hartung, J. S. and Levy, L. (2006). Quantitative real-time PCR for detection and quantification of *Candidatus Liberibacter* species associated with citrus huanglongbing. *Journal of Microbiological Methods*, **66**, 104–115.
- Loiseau, M., Renaudin, I., Cousseau-Suhard, P., Lucas, P. M. and Gentit, P. (2015). Assessment of performance criteria of *Candidatus Liberibacter solanacearum* detection methods. In: Testa-EPPO Conference on Diagnostics for Plant Pests. [https://www.eppo.int/MEETINGS/2015\\_meetings/wk\\_testa](https://www.eppo.int/MEETINGS/2015_meetings/wk_testa)
- Loiseau, M., Renaudin, I., Cousseau-Suhard, P., Lucas, P. M., Forveille, A. and Gentit, P. (2017a). Lack of evidence of vertical transmission of ‘*Candidatus Liberibacter solanacearum*’ by carrot seeds suggests that seed is not a major transmission pathway. *Plant Disease*, **101**, 2104–2109. <https://doi.org/10.1094/PDIS-04-17-0531-RE>
- Loiseau, M., Renaudin, I., Cousseau-Suhard, P., Poliakoff, F. and Gentit, P. (2017b). Transmission tests of ‘*Candidatus Liberibacter solanacearum*’ by carrot seeds. *Acta Horticulturae*, **1153**, 41–46.
- Mawassi, M. Dror, O., Bar-Joseph, M., Piasezky, A., Sjölund, J. M., Levitzky, N., Shoshana, N., Meslenin, L., Haviv, S., Porat, C., Katsir, L., Kontsedalov, S., Ghanim, M., Zelinger-Reichert, E., Arnsdorf, Y. M., Gera, A. and Bahar, O. (2018). ‘*Candidatus Liberibacter solanacearum*’ is Tightly Associated with Carrot Yellows Symptoms in Israel and Transmitted by the Prevalent Psyllid Vector *Bactericera trigonica*. *Phytopathology*, **108**, 1056–1066. <https://doi.org/10.1094/PHYTO-10-17-0348-R>
- Oishi M., Hoshino S., Fujiwara Y., Ushiku S., Kobayashi Y. and Namba, I. (2017). A comparison of protocols to detect *Candidatus Liberibacter solanacearum* from carrot seeds, research on the effectiveness of propidium monoazide treatment and evaluation of seed transmission in carrot seeds. *Research Bulletin of the Plant Protection Service Japan*, **53**, 111–117.
- Teresani, G. R., Bertolini, E., Alfaro-Fernández, A., Martínez, C., Tanaka, F. A. O., Kitajima, E. W., Roselló, M., Sanjuán, S., Ferrándiz, J. C., López, M. M., Cambra, M. and Font, M. I. (2014). Association of ‘*Candidatus Liberibacter solanacearum*’ with a vegetative disorder of celery in Spain and development of a real-time PCR method for its detection. *Phytopathology*, **104**, 804–811.